

Space requirements













(unit:mm)

Product specifications

	((i)) Cell³iMager	(♠) Cell³iMager neo				
Well plates	6 to 384 wells, flat and U-shape bottom types and 35 mm dish Microscope slide, 25cm² flask (only for Cell³iMager neo)					
Plate number	4	1				
Scan speed	54 seconds per plate (2400 dpi high speed mode)	50 seconds per plate (2400 dpi high speed mode)				
Automated measuring time	20 seconds per plate (2400 dpi high speed mode)					
Scan resolutions	200dpi, 300dpi, 600dpi, 1200dpi, 2400dpi, 4800dpi, 9600dpi					
Light source/optical system	White LEDs/transmission illumination					
Measuring	Number, area, estimated volume, contour, optical density, circularity and edge sharpness					
Output format	- Captured image display and numerical data confirmation using base software - Granular distribution and proliferation trend chart confirmation - Output of measured data in CSV file					
Environment	Recommended: Temperature 18 to 28 °C(64 to 82°F)					
Power requirements	Single phase 100 - 240V 0.12 kw 0.9A	Single phase 100 - 240V 0.115 kw 0.85A				
Dimensions (W x D x H)	492 x 755 x 294 mm	4 mm 594 x 381 x 294 mm				
Weight	Approx. 40 kg	Approx. 25 kg				

Specifications and equipment designs are subject to change without notice.

SCREEN Holdings Co., Ltd.

KYOTO(Head office) / Tenjinkita 1-1, Teranouchi-agaru 4-chome, Horikawa-dori, Kamigyo-ku, Kyoto 602-8585, Japan

Life Science Business Development and Sales Division

KYOTO(Rakusai) /

Furukawa-cho 322,Hazukashi,Fushimiku,Kyoto 612-8486, Japan Phone : +81-75-931-7824 / Fax : +81-75-931-7826

Te Yasukuni-Kudan-Minami Bldg., 2-3-14 Kudan-Minami Chiyoda-ku, Tokyo, 102-0074 Japan Phone: +81-3-3237-3950 / Fax: +81-3-3237-3938


www.screen.co.jp

BDS-102(E)(R2) Printed in Japan 12-15.010BB





High-Speed 3D Cell Scanner



A Bright Field Spheroid Counter, Ideal for Label-Free Assays of 3D Cell Culture

Compatible with various 3D culture plates, the Cell³iMager allows the use of 6-well to 384-well plates and a 35 mm dish.

It comes with an extensive and convenient software which enables the reading and image composition of soft agar medium, displaying both a full view of a well plate and magnified view of a well center, and the editing of measuring recipes.


A Compact Spheroid Counter with Greater Scanning Operability Without Compromising the Functions of Cell³iMager

Cell³iMager's sister unit, Cell³iMager neo has a compact design enabling the reduction of footprint by 40%. In addition to the existing functions of Cell³iMager, this highly versatile 3D spheroid counter provides excellent performance and scalability with even faster scanning speed for a single well plate and its numerous new software.









Perfect for a lab desk

Fearuves of high-speed 3D cell scanner

Point 1

High-speed well plate scanning and automatic measuring

Scanning 3D culture plates in just 54 seconds per plate, it enables high-speed measuring of the area, estimated volume, diameter and number of spheroids.

Point 2

Reagent-free and label-free easy imaging system

The label-free system, eliminating a series of processes including reagent dropping, stirring and fluorescence measuring, enables faster operation with fewer processes.

Point 3

High correlativity with drug susceptibility evaluation of ATP reagents and excellent spheroid image capturing performance

Installed software precisely captures spheroids by minimizing the effect of shadows created by meniscus on the peripheral area of wells. This excellent image capturing performance achieves high correlativity with drug susceptibility evaluation of ATP reagents.

Point 4

Follow-up and morphological observation of the same well allows time course observation and effortless image saving



Point 5

(Cell'iMager neo

Simple and user-friendly operability

Plates by all manufacturers and a wide selection of sample recipes of each well plate type are available.

Adopted graphical user interface (GUI) enables each operation step, based on well plate and scanning images.

Point 6

Designed to accommodate the automation of continuous scanning, it facilitates high-throughput observation

With the automated configuration of a platform, arm, dispensing heads, grip-hand and carrier, a sequence of processes; transferring well plates, dispensing fluid, scanning and storing is automated (under development). It is ideal for both time course observation and high-throughput scanning.



High-speed well plate scanning and automatic measuring

Scanning 3D culture plates in approximately 1 minute per plate, it then prepares a recipe for measuring based on the scanned images. In addition to our standard recipe samples, researchers can easily prepare a recipe that they want by fine adjustments of such sample recipe. Once the recipe is ready, measurement can begin, which takes approximately 20 seconds. As the follow-up stage involves scanning and auto-measuring only, an entire operation can be completed in approximately one and half a Follow-up run



0.5 min.

Scanning

Following completion of auto-measuring, it proceeds to scanned image observation and measured numerical data analysis by converting to graph form.















Spheroids can be observed in a raw image high quality TIFF format. Fast and easy image viewing is possible by selecting a whole well plate or individual well in a single image. The results obtained by auto-measuring are converted to graphs showing the proliferation trend and spheroid distribution, thereby aiding trend analysis of drug susceptibility.

All imaging photographs and test results can be backed up and restored by interconnecting. Various computed numerical data is output in the CSV format, enabling import into a spreadsheet or other user familiar software format.




Supporting the Reading and Image **Composition of Soft Agar Medium**

Multiple images can be obtained by changing the Z-direction.

2. A single composite image from obtained multiple images Two composite methods:

Prioritize pixel valuePrioritize edge contrast

Cells cultured in the soft agar medium are found in its overlapping layers. A new function allows the compositing of one image by imaging multiple layers in the Z-direction. Two modes; prioritize pixel value and prioritize edge contrast can be selected for computing the area and number of cell masses in the soft agar medium more precisely.




ATP assay

Reagent-free and label-free easy imaging system


Faster processing / Fewer processes

The optical label-free system enables easy operation by simply placing culture plates on the scanning table for scanning and measuring. Reagent dropping, stirring and subsequent fluorescence measuring are no longer necessary, allowing drastic time reduction and simplified operation.



High correlativity with drug susceptibility evaluation of ATP reagents and excellent spheroid image capturing performance

High correlativity has been confirmed between estimated volume measured by Cell3iMager and the relative light unit under the ATP assay, for drug susceptibility tests carried out on 50 cell types. This assures that it can be also used as a plate reader to test cell proliferation and activity.



Cell types of which the correlativity with ATP assay has been confirmed*

[Cell types]				
HT29	Colo201	Li-7	PP5	Capan-2
MiaPaCa-2	HCT116	SUIT-2	A549	WiDr
BxPC-3	RKO	JHH-5	SW620	DU145
T47-D	MCF-7	Panc-1	SKOV-3	NCI-H165
HLF	MM231	BT474	SW480	DLD-1
JHH-7	Capan-1	OVCAR-3	SK-BR-3	huH-1
HCC1954	HCC1806	BT20	HepG2	
HeLa	MKN-74	HLE	MKN-45	
AsPC-1	Colo205	ACHN		
JHH-2	NCI-H1993	3T3L1		
NCI-H441	NCI-H23	Hep3B		
Hs578T	NCI-H1975	PC-3		


*By inhouse tests

In order to obtain high correlativity with ATP reagents, software is installed which captures spheroids precisely. Cell³iMager uses the transmission illumination system for estimation of 3D cultured spheroid volume, which may produce the meniscus effect on the image depending on the culture medium. The software installed minimizes such shadow effect created by menisci and helps capture spheroids at the peripheral area of wells,

Compatible with various 3D culture plates


on their background color intensity, preventing image capturing performance

failure resulting from the shadow.







Follow-up and morphological observation of the same well allows time course observation and effortless image saving.



Time course statistics

Different culture environments within a well plate and compound test patterns can be set per group. This allows not only the tracking of changes and development of images based on data imaged over time, but also the numerical control of changes in the number, area and volume of cells. Further, with the function that converts data to graph format by the single touch of a button, the researcher can have a full view of his/her desired data e.g. the timing of drug effects.






▲The plate map shows	changes/development by icons.
----------------------	-------------------------------


	Name	2days(48)	3days(72	4days(9€	6days(1	9days(2	11days	13days	Comment1	Comment2
V	ATEL	552053	66427	77709	72112	10271	1152	1306	LNCap / PRMI	1000nM
V	C1:D4	513516	65549	79966	14867	30840	3267	3356	LNCap / PRMI	100nM
V	E1F4	505887	70498	11248	27193	39782	3772	3464	LNCap / PRMI	10nM
V	G1:H4	510535	65986	11603	27222	36315	3359	3104	LNCap / PRMI	Control
	A5 88	997666	91791	10519	12167	11308	1292	1491	MCF7 / RPMI1	1000nM
	0508	107829	89095	51278	13142	20159	2093	1877	MCF7 / RPMI1	100nM
	55.75	109317	11174	15537	23183	21200	1975	2179	MCF7 / RPMI1	10nM
	G5 H8	104372	11516	14855	21374	18990	1854	2190	MCF7 / RPMI1	Control
	49:812	139953	10548	85369	10735	36663	4081	4273	3T3 rasLacZ/	1000nM
	C9:D12	160232	12513	11508	20702	36340	4532	5804	3T3 rasLacZ/	100nM
	E9:F12	162166	15696	19553	35566	48903	6029	7168	3T3 rasLacZ/	10nM
	G9:H12	155781	17564	23032	35657	47967	5731	6726	3T3 rasLacZ/	Control

▲Time course data sheet of the spheroid area by group (Inclusion of comments re. cell types and culture medium volume is possible.)

Selectable scanning resolution

Select 2400 dpi resolution to obtain rough drug susceptibility at high speed, and 4800 dpi or 9600 dpi for more precise observation as needed





*It is possible to scan the entire well plate or a specific well.


- 2400 dpi: Approximately 10 µm
- 4800 dpi: Approximately 5 µm
- 9600 dpi: Approximately 2.6 µm


Point

Simple and user-friendly operability



Set the size, concentration, circularity, etc. and select the spheroids subject to measuring.





Granular distribution of specified wells, a graph of the follow-up proliferation development and other test data can be created using the standard software.

Detail analysis on computer

same folder with measuring

recipes and measured data for easy access as needed.

.


Storage

Measured data output in the CSV format can be used for detail analysis on a spreadsheet or other software format. The spheroid area, number and estimated volume are computed for use in CSV data and graphs.

Fibriform debriis separation function



A new function is added which enables the measurement of cell spheroids even when the background density differs, by defining the volume of the peripheral area of wells. It separates subject spheroids from overlapping fibriform debris and improves the precision of measurement. As both the full view of a well plate and magnified view of a single spheroid's center are available, work can be executed with a comprehensive view of the subject.







Through collaboration with a liquid handling robot manufacturer, a sequence of processes; transferring well plates from a carrier, dispensing fluid, scanning with the spheroid counter and storing by a carrier is automated (under development), with the configuration of a platform, arm, dispensing heads, grip-hand and carrier. It enables mass continuous processing, while providing optimal efficiency for both time course observation and high-throughput scanning.

